


IMATINIB IMPURITY E

<p>Chemical name: Imatinib dimer; Des(methylpiperazinyl-N-methyl) Imatinib Dimer Impurity; 1,4-Bis-[4-[4-methyl-3-[[4-(pyridin-3-yl)pyrimidin-2-yl]amino]phenyl] carbamoyl] benzylpiperazine</p>	<p>CAS No. : [1365802-18-1]</p>
<p>Category : Impurities / Metabolite</p>	
<p>Molecular Formula : C₅₂H₄₈N₁₂O₂</p>	<p>Structure:</p> 
<p>Molecular Weight : 873.02</p>	
<p>Product Code : RBPL-I0101</p>	