


RABEPRAZOLE IMPURITY A

<p>Chemical name: Rabeprazole sulphone; 2-[[[4-(3-Methoxypropoxy)-3-methyl-2-pyridinyl]methyl]sulfonyl]-1H-benzimidazole; Rabeprazole Related Compound D</p>	<p>CAS No. : [117976-47-3]</p>
<p>Category : Impurities / Metabolite</p>	
<p>Molecular Formula : C₁₈H₂₁N₃O₄S</p>	<p>Structure:</p> 
<p>Molecular Weight : 375.44</p>	
<p>Product Code : RBPL-R0202</p>	